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### NÃO INJECTAR.

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## Instructions for Use

### DO NOT INJECT.

Vitigel™ Surgical Hemostat ("Vitigel") must not be injected into blood vessels.

### DEVICE DESCRIPTION AND ACTION

Vitigel includes:

1. A Vitigel syringe containing a suspension of 20 mg/mL bovine collagen and at least 300 U/mL bovine thrombin (THROMBIN-JMI®) in a 40 mM CaCl<sub>2</sub> buffer.
2. A Delivery System containing accessories used for joining the Vitigel syringe to the autologous plasma syringe, and for mixing and delivering the fluids.
3. A Transfer Syringe used to receive the patient's plasma and provide a means for placing the plasma into the sterile operating field.

Vitigel is used in conjunction with the patient's own plasma. A product called CellPaker® Collection Device is supplied separately and is used for obtaining the patient's plasma.

The bovine thrombin in Vitigel converts autologous plasma fibrinogen to fibrin, which, in the presence of collagen, forms a collagen/fibrin gel matrix which adheres to the bleeding site. Vitigel is biocompatible, non-pyrogenic, and resorbed within 30 days and is intended to be left *in situ*.

**INDICATIONS:** Vitigel is indicated in surgical procedures (other than in neurosurgical and ophthalmic) as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

### CONTRAINDICATIONS

- Do not inject Vitigel into blood vessels or allow it to enter blood vessels. Extensive intravascular clotting and even death may result.
- Do not use Vitigel in patients known to be sensitive to any of its components and/or to materials of bovine origin.

### WARNINGS

- Vitigel is not intended as a substitute for meticulous surgical technique and the proper application of ligatures or other conventional procedures for hemostasis.
- Vitigel should not be used in the presence of infection. Vitigel should be used with caution in contaminated areas of the body. If signs of infection or abscess develop where Vitigel has been applied, reoperation may be necessary in order to remove the infected material and allow drainage.
- Excess Vitigel may be removed by gentle irrigation from the site of application when used in, around, or in proximity to foramina in bone, areas of bony confine, the spinal cord, and/or the optic nerve and chiasm.
- The safety and effectiveness of Vitigel for use in neurosurgical and ophthalmic procedures have not been established.
- Do not use Vitigel in instances of profuse arterial bleeding. It should not be used where blood or other fluids have pooled or in cases where the point of hemorrhage is submerged.
- Vitigel should not be used for controlling post-partum bleeding or menorrhagia.
- The safety and effectiveness of Vitigel have not been established in children and pregnant women.

### WARNING

The use of topical bovine thrombin preparations has occasionally been associated with abnormalities in hemostasis ranging from asymptomatic alterations in laboratory determinations, such as prothrombin time (PT) and partial thromboplastin time (PTT), to severe bleeding or thrombosis which rarely have been fatal. These hemostatic effects appear to be related to the formation of antibodies against bovine thrombin and/or factor V which in some cases may cross react with human factor V, potentially resulting in factor V deficiency. Repeated clinical applications of topical bovine thrombin increase the likelihood that antibodies against thrombin and/or factor V may be formed. Consultation with an expert in coagulation disorders is recommended if a patient exhibits abnormal coagulation laboratory values, abnormal bleeding, or abnormal thrombosis following use of topical thrombin. Any interventions should consider the immunologic basis of this condition. Patients with antibodies to bovine thrombin preparations should not be reexposed to these products.

### PRECAUTIONS

- Vitigel is supplied as a sterile product for single use only. Do not resterilize.
- Once hemostasis is achieved, any excess Vitigel may be carefully removed by gentle irrigation.
- As with other hemostatic agents, do not aspirate Vitigel into extracorporeal cardiopulmonary bypass circuits or autologous blood salvage circuits. It has been demonstrated that fragments of collagen based hemostatic agents may pass through 40µ transfusion filters of blood scavenging systems.
- Vitigel should not be used in conjunction with methylmethacrylate or other acrylic adhesives. Microfibrillar collagen has been reported to reduce the strength of methylmethacrylate adhesives used to attach prosthetic devices to bone surfaces.
- Vitigel should not be used for the primary treatment of coagulation disorders.
- The safety and effectiveness of the combined use of Vitigel with antibiotic solutions or powders have not been established.
- The safety and effectiveness of Vitigel for use in urological procedures have not been established through randomized clinical study.
- In urological procedures, Vitigel should not be left

in the renal pelvis or ureters to eliminate the potential foci for calculus formation.

Pregnancy - Category C - Animal reproduction studies have not been conducted with THROMBIN-JMI. It is also not known whether THROMBIN-JMI can cause fetal harm when administered to a pregnant woman or can affect reproduction capacity. THROMBIN-JMI should be given to a pregnant woman only if clearly indicated.

**ADVERSE EVENTS:** In a randomized prospective, concurrently controlled clinical trial, 318 patients were treated with Vitigel or the Control (a collagen absorbable hemostat for the non-cardiac procedures, and other conventional methods selected by the surgeon for cardiac procedures). The most common adverse events recorded during and after the application of the hemostatic agents were fever, pain, nausea, and atelectasis. None of the adverse events were judged to be attributed to Vitigel.

### DIRECTIONS FOR USE

#### HOW TO ASSEMBLE VITAGEL SURGICAL HEMOSTAT

Refer to Figure 1 and use the following steps to assemble Vitigel Surgical Hemostat.

- 1.1 Remove CellPaker from the centrifuge and hold in a vertical position with the Luer tip up. Remove cap from Luer tip. Turn up the red outersleeve of CellPaker slowly until air is expelled.
- 1.2 Connect CellPaker to Transfer Syringe. Remove and discard transfer syringe Luer cap. Transfer syringe tip is sterile. Attach transfer syringe to CellPaker.
- 1.3 Transfer Plasma from CellPaker to Transfer Syringe. Re-attach the plunger to the CellPaker by turning clockwise and use it to transfer plasma into the syringe. The volume of plasma transferred should be equal to the volume of the Vitigel syringe being used.
- 1.4 Place Plasma Containing Transfer Syringe into Sterile Field. The transfer syringe containing plasma is sterile. Circulating nurse peels open transfer syringe pouch. Scrub nurse unscrews and removes plasma syringe and places into sterile field.
- 1.5 Place the Sterile Vitigel Syringe into the Sterile Field. The Vitigel suspension syringe is sterile. Circulating nurse peels open the Vitigel syringe pouch. Scrub nurse unscrews and removes the Vitigel syringe and places into the sterile field.
- 1.6 Expel Air from Vitigel and Plasma Syringes.
- 1.7 Assemble Syringes with Delivery System. It is important to attach the cannula or sprayhead last.
  - a. Attach Luer ends of Vitigel and plasma syringes to joiner.
  - b. Align ends of plunger rods until even. Expel syringe contents, if necessary.
  - c. Slide syringe clip over the ends of plunger rods.
  - d. Slide support over assembled device until firmly seated.
- 1.8 Attach Spray Head, Cannula or an equivalent delivery device cleared by the FDA for use with Vitigel according to the surgeon's preference. Specific instructions for the application of Vitigel with devices other than the Spray Head or Cannula are provided with each cleared device.

The Vitigel Surgical Hemostat is now ready for use.

#### HOW TO APPLY VITAGEL TO THE BLEEDING SITE

Refer to Figure 2 and use the following steps to apply Vitigel to the bleeding site.

- 2.1 Discontinue all blood recovery and cell-saving devices, as appropriate. Blot surface dry.
- 2.2 Apply product as a thin, uniform coating and overlap the edges to ensure complete coverage. *Avoid disrupting the gel.*
- 2.3 If discrete bleeding is observed, additional Vitigel may be applied by underlying the formed gel with the tip of the cannula and delivering the material directly to the site.  
If bleeding continues, remove the first application as completely as possible, blot, and reapply Vitigel.

**HOW SUPPLIED:** Vitigel Surgical Hemostat is supplied as a single use only 5 mL treatment syringe, which is packaged with sterile delivery components and a sterile transfer syringe. Vitigel is combined with an equal volume of patient plasma.

Sterile Delivery System and sterile Transfer Syringe are included.

*Do not use if package is damaged or if syringe caps are dislodged or missing.*

The CellPaker Collection Device is available separately in multiple unit boxes.

Vitigel and its accessories contain no natural rubber latex.

**STORAGE CONDITIONS:** Store Vitigel Surgical Hemostat at 2 - 8 °C. **Do not freeze.**

**CAUTION: RX USE ONLY – Limited by Federal law to use by or on the order of a physician.**

**U.S. PATENT NUMBER:** Vitigel Surgical Hemostat may be covered by one or more of the following: U.S. Patent No. 5,290,552; U.S. Patent No. 6,096,309; and 5,997,811.

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© 2007, Orthovita®, Inc.